## User guide

RMJ Health AB  
Telephone: +46 708-51 72 61  
E-mail: [kontakt@rmjhealth.com](mailto:kontakt@rmjhealth.com)

**Index**

1. Warnings .....2

2. Intended use .....5

2.1 Before use .....5

2.2 After use .....6

3. Example of exercises .....6

3.1 Biceps.....7

3.2 Triceps.....8

3.3 Shoulders.....9

3.5 Back .....10

3.6 Chest.....11

4. Warranty and lifetime .....11

## 1. Warnings

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### Thank you for purchasing a product from RMJ Health AB!

Our focus is on high quality of training and rehabilitation. In our product range you will find everything from knee pads to exercise bands.

**Warning:** Before you begin your training or rehabilitation, it is important to first read this. Always consult a physical therapist or personal trainer if you are unsure how an exercise should be performed. If you experience pain or feel faint during exercise, stop training immediately and consult a physician. The product is suitable for home use as well as at a sport- or training facility. Always inspect the product and its point of attachment before use. You must always be absolutely sure that the anchor point and the product are stable enough to withstand the load before performing exercises. For questions, contact our support via mail at, [support@rmjhealth.com](mailto:support@rmjhealth.com)

**Warning:** This product is not a medical article, but should be used as an adjunct to facilitate physical activity. If you have a medical condition you should contact a doctor before use. RMJ Health Ltd cannot guarantee that damage will not occur even if the product is used correctly. With proper use, the right equipment, common sense and knowledge of the risks the risk will be minimized.

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## 2. Intended use

Activation Grip is a glove, to which you can attach a weight to be able to train the upper body strength without holding an object in your hand. Activation Grip can also be used to train the strength of the upper extremity by imitating the sport you practice as much as possible.

Intended use for Activation Grip is as an aid in the rehabilitation of people who can't hold an object in his/her hand such as hand injured individuals and individuals with impaired function of the hand as rheumatism. The product can also be used by athletes and exercisers who usually hold their sports gear in the hand.

The use of Activation Grip is intended in the home, fitness center or sports center. Follow the instructions carefully before using the product in order to minimize the risk of accidents.

Maximum load for Activation Grip is 45 kg.

### 2.1 Before use

Follow the bulleted list each time before using the Activation Grip:

- Make sure the glove is not damaged by following these points:
  - Is the eyelets, rings and cord damaged?
  - Is the resistance connected properly in the glove?
  - Is the resistance point of attachment stable and stays in place for training?
  - Is the glove correctly placed on your hand? The print should be on the back of the hand. (See image below)

Figure 1: Suggested coupling Activation Grip to a cable pulley machine.

## 2.2 After use

Follow the bulleted list each time after using the Activation Grip:

- Make sure the glove is not damaged by following these points:
  - Is the eyelets, rings and cord damaged?
  - If you share the glove with other people, make sure to clean it with fragrance free soap and water after use to minimize the risk of infection.

## 3. Example of exercises

*Activation Grip at the gym:*

Connect the Activation Grip to a cable pulley machine by putting the snap hook in the ring on the inside / outside of the hand, depending on the exercise. Select the desired resistance of the machine.

*Activation Grip at home:*

Connect Activation Grip to the exercise band. Choose a secure point of attachment in the home that cannot tip over or break. Door or banisters are recommended. Doorstop is available as an accessory on [www.activationgrip.se](http://www.activationgrip.se)

Below are some simple exercises that provide a good workout for the upper body. For more exercise suggestions, visit [www.activationgrip.se](http://www.activationgrip.se)

### 3.1 Biceps

*At the gym:*

*Start*

*Stop*

Move the hand against the shoulder as far up as you can. Proceed slowly back to the starting position and repeat the movement.

*At home:*

*Start*

*Stop*

Move the hand against the shoulder as far up as you can. Proceed slowly back to the starting position and repeat the movement.

### 3.2 Triceps

*At the gym:*

*Start*

*Stop*

Keep your arm in a 90 degrees angle and extend your arm completely. Proceed slowly back to starting position and repeat.

*At home:*

*Start*

*Stop*

Hold up your arm with your elbow pointing forward. Stretch out your arm as straight as possible while maintaining the same angle of the upper arm. Proceed slowly back to starting position and repeat.

### 3.3 Shoulders

*At the gym:*

*Start*

*Stop*

Start with the arm along the side of the body and then lift with your arm stretched straight out from the body at shoulder height. Slowly return to start position and repeat. NOTE: If you feel any the pain in the shoulder area, do not go higher then what is painless!

*At home:*

*Start*

*Stop*

Start with the arm along the side of the body and then lift with your arm stretched straight out from the body at shoulder height. Slowly return to start position and repeat. NOTE: If you feel any the pain in the shoulder area, do not go higher then what is painless!

### 3.5 Back

*At the gym:*

*Start*

*Stop*

Stand with your legs slightly bent and arch your back. Pretend you are holding a stick and pull your hands up towards your chest. Slowly return to starting position and repeat.

*At home:*

*Start*

*Stop*

Stand slightly leaning forward with your weight on your front leg. Pull your arm back, as high as possible. Remember to also use the back muscles in the movement. Slowly return to starting position and repeat.

### 3.6 Chest

*At the gym:*


*Start*


*Stop*

Start in a cross-like position where the arms are held straight out from the body. Push your hands together in front of your body still with outstretched arms. Remember to tighten up the chest. Slowly return to starting position and repeat.

*At home:*

*Start*

*Stop*

Start with hands parallel with your shoulders. Stretch your arms straight out. Remember to tighten up the chest. Slowly return to starting position and repeat.

## 4. Warranty and lifetime

RMJ Health offer 1 year warranty from date of purchase on Activation Grip. The warranty covers the production & manufacturing of the product and does not cover damage to the product due to improper use by the user.

Activation Grip has an estimated lifetime of 1 year in regular exercise.